

Ambu[®] Man Basic



Key features

- Cost efficient BLS training manikin
- Robust and reliable
- Realistic training of chest compression
- Visualisation of the compression depth
- Ambu hygienic system für mouth-to-mouth and mouth-to-nose ventilation
- Can be used for bag-mask ventilation
- Placement of Guedel tubes (When using the Guedel face piece)
- Realistic ventilation resistance
- Visualisation of the correct ventilation volume
- Visible chest rise during ventilation
- Anatomically correct closed Torso
- Easy AED training

AmbuMan Basic

The AmbuMan Basic is a cost efficient training manikin that is designed to fulfill all requirements for the BLS algorithm training.

The closed shaped torso with its anatomically correct landmarks trainees can also learn and train the use of Automatic External Defibrillators (AED) and the placement of pads.

The AmbuMan Basic allows an optimal training of chest compression and ventilation with realistic resistance and expiration through mouth and nose. The unique hygienic system prevents cross contamination. The correct compression depth and required ventilation volume can be monitored by the feedback indicator on the manikin.

The AmbuMan Basic is also available as 4 pack to enable group training.

Realistic ventilation

The students will feel the realistic airway resistance and the chest will start rising when the head is correct tilted and air blown into the manikin.

Being able to look, listen and feel according to BLS guidelines.

Kosteneffizient and reliable

The Ambu Man Basic is made to be cost efficient and reliable for the BLS training.

Class training

For an easy transportation & set up possibility the AmbuMan Basic comes in a soft bag either packed as single manikin or in set of four manikins.

Specifications

Weight Torso	6 kg
Length	70 cm
Width	37 cm
Weight 4 x Torso and Bag	30 kg
Max. Ventilation volume	1200 ml
Max. compression depth	7,5 cm

Materials

Basic unit	ABS (Acrylnitril-Butadien-Styrol)
Head	Polyvinylchlorid, hart
Chest skin	Polyvinylchlorid, weich
Face piece	Polyvinylchlorid, weich
Headbag, disposable	Polyethylen
Trainingsuit	Polyester
Carrying bag	PVC-beschichtetes Nylon

Order information

Item no	Description
A268 005 000	AmbuMan Basic inkl. training jacket, carry case, 1 face piece and 25 head bags
A268 006 000	AmbuMan Basic 4 pack inkl. 4 AMBuMAN Basic Torso with training jacket, carry case, 4 face pieces and 100 head bags
Accessories	
234 000 702	Ambu head bags, disposable (set of 100)
234 000 703	Ambu face piece (Set of 5)
234 000 725	Ambu Guedel face piece



AmbuMan Basic



Hygienic System



Compression indicator

Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com