

## DATASHEET


---

Ambu® Subdermal Hooked Needle and Ambu® Subdermal Curved Needle are stainless steel electrodes, which are designed to stay in place during usage. The electrodes have color coded lead wires which aid the identification of electrodes. The products come with a 1.5 mm touch proof DIN 42802 connectors and can be used with most multi-channel nerve monitors.

The hooked and curved needle electrodes are developed to be used for monitoring, recording and stimulation of bio- and nerve potential signals.

# Ambu® Subdermal Hooked Needle

# Ambu® Subdermal Curved Needle



## KEY FEATURES

Designed to stay better in place during usage

Color coded lead wires aid the identification of electrodes

Needle electrode made of stainless steel

1.5 mm touch proof DIN 42802 connector

---

## RECOMMENDED APPLICATIONS

Intraoperative Monitoring (IOM)

Electroencephalogram (EEG)

Electromyography (EMG)

Evoked Potentials (EP)

# Ambu

## SPECIFICATIONS

### Electrode

Connector type	1.5 mm touch proof DIN 42802
Sterilization method	Ethylene oxide (EO)

### Environment

Electrode, lead wire and packaging are not made with natural rubber latex
Packaging is PVC free


## MATERIALS

### Electrode


Needle	Stainless Steel
Lead wire	PVC insulated Tin Plated Copper Wire
Shrink tube	Polyolefin
Protection tube	Silicone

### Packaging

Pouches, transparent layer	Polyester/Polyethylene
Pouches, paper layer	Paper
Boxes	Cardboard



Subdermal Hooked Needle



Subdermal Curved Needle

## ORDERING SPECIFICATIONS

Item no.	Lead wire length		Lead wire color	Shelf life in months (unopened pouches)	Units/pouch	Packaging	
	cm	inch				Units/inner box	Units/outer box
<b>Subdermal Hooked Needle</b>							
74912-120/24	120	47"	○ ● ● ● ● ●	36	1	24	1152
74912-200/24	200	79"		36	1	24	1152
74912-250/24	250	98"		36	1	24	1152
<b>Subdermal Curved Needle</b>							
75021-150/25	150	59"	○ ● ● ● ●	36	1	25	1200
75021-200/25	200	79"		36	1	25	1200

# Ambu

#### Distributed by:

**Ambu A/S**

Baltorpbakken 13  
2750 Ballerup  
Denmark  
T +45 72 25 20 00  
ambu.com

#### Manufactured by:

**Technomed Europe**

Amerikalaan 71  
6199 AE Maastricht-Airport  
The Netherlands

CE<sub>0344</sub> US: Rx only